

NDC Translator™

NDC Translator™ software translates the raw data from an NDC barcode into an NDC number format usable by applications such as BPOC, BCMA and NDC billing. NDC Translator™ is compatible with hospital information systems and most models of scanners.

How It Works:

After scanning the bar code, the NDC Translator reviews the raw data from the bar code scanner. If it finds an NDC number it sends it to the application in one of four user selected output formats:

- 1) The 10 digit format of the NDC number
e.g. 0045064065
- 2) The 11 digit (zero added) format of the NDC number
e.g. 00045064065
- 3) The 10 digit format of the NDC number with dashes
e.g. 0045-0640-05
- 4) The 11 digit (zero added) format of the NDC number with dashes
e.g. 00045-0640-65



This provides the user with quick entry of the NDC number needed for accurate billing and documentation of the product used.

It also can be set to output lot number / expiration date information and when the expiration is included in the bar code, **it will pop up an expired warning.**

Who:

Anyone who needs to input NDC numbers into their software applications, such as for:

- BPOC / BCMA
- Dispensing Documentation
- NDC Billing Compliance
- Drug Repackaging
- Drug Return and Disposal
- Inventory Management

www.RxScan.com

Sales@RxScan.com • 1-800-572-2648

©2009 RxScan, Ltd. All Rights Reserved.

RxScan

NDC Translator™